



 **smith&nephew**
REDAPT[◇]
Revision Femoral System











The power to *redapt*
any revision





REDAPT[®] Revision Femoral System by Smith & Nephew is uniquely designed to address the challenges in today's revision hip arthroplasty

Paprosky Femoral Revision Classification

Type 1	Type 2	Type 3A	Type 3B	Type 4
				
				

One system for every revision

There's nothing routine about revisions. REDAPT allows you to address every patient effectively and efficiently.

ROCKTITE[®] fixation

Patented multi-level flutes and a three-degree taper design provide rock-solid distal fixation in all bone types.



Stem cross section

- Fixation in various bone types
- Independent achievement of joint stability
- Proximal/distal mismatch without compromising stem strength
- Predictable stem position
- Surgical efficiency

Adaptable implant design

The adaptable design allows you to customize the proximal/distal mismatch for the best stem fit without compromising stem strength.

Further, you are able to achieve joint stability through independent adjustment of neck angle, height, offset and version.



Reproducible stem position

The REDAPT[®] system's reproducible stem position is achieved by the combination of ROCKTITE[®] fixation and efficient, easy to use instrumentation.



Smith & Nephew, Inc.
7135 Goodlett Farms Parkway
Cordova, TN 38016
USA

www.smith-nephew.com

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders and Inquiries: 1-800-238-7538

°Trademark of Smith & Nephew

©2012 Smith & Nephew, Inc.
All rights reserved.
7138-1685 REVA 10/12